

# HuMap<sup>™</sup> + Neurotransmitter (NT) testing

- Patients are complex. Similar symptoms can be the result of many different, and individualized, neurotransmitter and hormone imbalances.
- Mental health problems are not isolated in the brain alone. There are many factors affecting our mental wellbeing and by assessing neurotransmitters and reproductive and adrenal hormones we can get a more well-rounded understanding of how nutrition, digestion, inflammation, stress, metabolism, and methylation may be contributing to the symptom picture.
- Test in order to determine which therapeutics are appropriate.





# HuMap™

### Hormone & Urinary Metabolites Assessment Profile


- Extensive estrogens and metabolites: E1, E2, E3, 2-OH-E1/E2, 4-OH-E1/E2, 16-OH-E1, 2-Methoxy-E1/E2,4-Methoxy-E1/E2
- Comprehensive metabolites of progesterone and androgens
- Diurnal cortisol and cortisone
- Metabolized cortisol, cortisone, corticosterone
- Key enzyme activity: COMT-11 $\beta$ HSD2, aromatase, 5 $\alpha$ -reductase
- Biomarker for oxidative DNA damage: 8-OHdG
- Abbreviated profiles available to monitor treatment response



Ne	eur	o Ba	isic	or	Compre	hensiveNeu	iro	tran	sm	hitter		
		_	_									
Neuro Basic Profile	; urine				DECTORS DATA	Comprehensive Neu	rotransr	nitter; urine			DCTORS DATA	
Order: SAMPLE REPORT Test: U12446-7800 Glient #: 12445 Doctor: Samjei Doctor Doctor's Data, Inc. 3755 Illinois Ave. St. Chartes, IL 60174	Age: 6 Sex: Fr		:20.0	Sample Collect Date Collect Wake Up Ti Collection I Date Receiv Date Report	me 07:30 Period 1st moming void ved 03/03/2022	Order: SAMPLE REPORT Test: U12365-7890 Citert #: 12345 Doctor: Sample Doctor Doctor's Data, Inc. 3755 limota Ave. St. Charles, IL: 60174	Age: Sex: Body	nt: Sample Patient 10 Female Mass Index (BMI) Mass Index (BMI)		Sample Collect Date Collect Wake Up Ti Collection I Date Recein Date Report	me 07:30 veried 1st morning void ed 03/03/2022	
Analyte	Result Un	nit per Creatinine	L	WRI	H Reference Interval			Init per Creatinine	1.1	WRI	H Reference Interval	
Serotonin	63.6	6/6H	4		60-125	Analyte Phenethylamine (PEA)	Result U		L	WRI	H Reference Interval	
Dopamine	120	hð <sub>i</sub> ð			125-250			nmol/g	-		42-160	
Norepinephrine	16.6	B,64	-		22-50	Tyrosine	48	µmol/g	-			
Epinephrine	1.5	hð <sub>i</sub> ð			1.6-8.3	Tyramine	3.4	µmol/g	- 4		2.8-8.5	
Norepinephrine / Epinephrine ratio	11.1		_	A	<13	Dopamine	258	6,64	_	A.	175-600	
Glutamate	14	nmol/g			12.0-45.0	3,4-Dihydroxyphenylacetic acid (DOPAC)	1520	hðyð		Δ	540 - 1850	
Gamma-aminobutyrate (GABA)	3.5	nmol/g		A	2.0-5.6	3-Methoxytyramine (3-MT)	151	nmol/g	_	*	122 - 278	
Glycine	694	nmol/g			450-2200	Norepinephrine	17.6	haya	A		29-69	
Histamine	14	<i>b</i> /64	Δ		14-44	Normetanephrine	195	PB/8		A	112-400	
Phenethylamine (PEA)	107	nmol/g			32-84	Epinephrine	1.9	hð <sub>i</sub> ð		-	2.1-14.5	
Creatinine	118	mg/dL			30-225	Metanephrine	84	haya		A	60 - 158	
Neurotransmitter Comments:						Norepinephrine / Epinephrine ratio	9.3		-	A	<13	
Urinary neurotransmitter levels provide an	overall assesse	ment of the body's a	hilly to make a	and break down never	transmitters and are	Tryptamine	0.4	µmol/g			0.65-1.6	
representative of whole body levels. Neuro systems. The enzymes, cofactors and pre-	stransmitters are	e secreted all through	the body in	neurons of both the o	entral and peripheral nervous	Serotonin	84.4	hð <sub>i</sub> ð	$ \Delta $		79-235	
nervous system. Therefore, alterations in a associated with many symptoms including	urinary neurotrai	insmitter levels asse	ssed in urine p	rovide important clini-	cal information, and may be	5-Hydroxyindoleacetic acid (5-HIAA)	5913	hð <sub>i</sub> ð	_	*	2500-9000	
pain. • Low range serotonin may contribute to mo						Glutamate	23	nmol/g		*	18.0 - 70.0	
serotonin may also be associated with poo lassitude. Failure to recenerate tetrahydrol	x sleep quality a	and appette change	s, as well as o	hronic fatique, meum	atoid arthritis, and over-all	Gamma-aminobutyrate (GABA)	5	nmol/g		Δ.	2.6-8.0	
lassitude. Failure to regenerate tetrahydrol could be reflected in urine. BH4 regenerati serotonin requires vitamin D, iron and vitar	ion may be supp	ported by folates, vit	tamin B3, C, m	olybdenum and zinc.	Additionally, production of	Glycine	966	nmol/g		A	700 - 2500	
theanine may affect serotonin function.						Histamine	12	hðyð			14-51	
<ul> <li>Low dopamine may be associated with ani increased addiction and other stimulation r</li> </ul>	seeking activitie	s. Failure to regene	rate tetrahydro	biopterin [BH4], an es	sential cofactor for departine	Taurine	428	µmol/g	4		420 - 1400	
synthesis, may decrease dopamine levels, molybdenum and zinc. Additionally, produc	and could be re tion of dopamin	affected in urine. BH ne requires vitamin <sup>1</sup>	44 regeneration D, iron and vita	may be supported b min B6. L-tyrosine, L-	y folates, vitamin B3, C, theanine and Mucuna pruriens	Creatinine	73.1	mg/dL		*	25-180	
may influence dopamine signaling. Low norepinephrine and low epinephrine magnetic section of the section of	nav be associat	ted with depression	and mood cha	nges as well as fatiou	e. difficulty concentrating.							
decreased ability to stary focused on tasks requiring vision (C. coper and nation (R3) - Upper range NIE ratio is consistent with po- h-methytranetasae (PMIT) request that interpretation in contox of cortical levelsh. Elevated phenethytamine (PRA) nay cont fermating at microbox. PRA and other to associated weth higher cortical levels. • Note: The reported low is low range mono quantities of required calactors. Further til	and diminished ). L-tyrosine, L-to or convention of requires SAMe, IPA axis function ribute to anxiety, e oxidase inhibit ace amines are	I sense of personally theanine and Mucur of norepinephrine to , magnesium and co n, with subsequent c , with very high leve fors or antipsychotic found in fermented	professional di paprofessional di epinephrine. ' artisol (adequa optimization of its having amp o medications, foods (wine, o	Ive. Norepinephrine is aence this pathway. This conversion is drive the HPA axis function where HPA axis function where heatamine-like effects. high protein diets, and heese, chocolate, etc.	: converted from dopamine en by the phenylethanolamine as cofactors. Suggest en clinically warranted. Elevations in PEA may occur a production by protein- .). Elevated PEA levels may be	Neurotransmitter Comments:     Irinary neurotransmitter Sevels provide an or system. The enymes, colladors and proceedings neurous system. Therefore, alterations in un associated with many symptoms including or pain. The enymes of the provide and the pro- sent of the procursor and the PEA, and via the procursor and so of the A. and via	ransmitters a insors in neuro inary neurotr ognitive and ted with depr	re secreted all throug rotransmitter metabol ansmitter levels asse mood concerns, dimi ression, attention defi	h the body, sm in gener used in urine hished drive cits and hyp	in neurons of both the o al are the same in the p a provide important clinic , fatigue and sleep diffic eractivity (ADHD), and t	ntral and peripheral nervous eriphery and in the central al information, and may be attes, cravings, addictions and loolar disorder. Phenvlalanine	DECTORS DATA

# Road Map

- Introduction to the HuMap™
- What are urinary hormone and metabolites?
- What is involved in metabolism of sex hormones?
- Major hormones and metabolites to pay attention to in your practice
- 8-OhdG and its role in oxidative stress
- Influences of hormones on neurotransmitters and vice versa the benefits of testing together



STE	EP 3: Collect Sample			
	ORE YOU BEGIN COLLECTING SAMPLE: te the Patient's Name, Date of Birth and		h Specimen Vial.	
	Day 1 Collection		Day 2 Collection	
	#1 Dinnertime: collect within 1-2 hours prior to eating	Midsleep: collect only	#3 Waking: collect within 10 mins of waking	
	#2 Bedtime: collect within 1 hour prior to sleep	if you wake to urinate while sleeping	#4 Post-waking: collect after 2-4 hours	
	Collect a midstream specimen of urine us If needed, you may reuse these cups. Pipette or pour urine into the properly lat	belled tube to the "fill line" mar	ked on the tube label.	
	Any remaining urine in the collection cup DO NOT urinate directly into the tube, wh			
	Close tube and rock gently.			
4	Record your NAME and the DATE and TIM	E OF COLLECTION on the tube	e, and on the Requisition Form.	
	Place the urine tube back into the origina Store the urine in the freezer with the froz		ssorbent pad.	7 DETORS SMALL SCIENCE - INSIGHT

				DISCONTINUE THE FOLLOWING
			Ð	Cortisol/glucocorticoid supplementation: - Certain medications such as asthma inhalers and hydrocortisone contain cortisol. To evaluate natural cortisol production, it is recommended to stop cortisol containing products 2 days prior to sample collection.
			5 days before and during test	Consult your practitioner for specific instructions before stopping any medications. Never discontinue prescription medications without first consulting your doctor.
	AT DAYS OF THE MONTH SHOULD I Women with menstrual cycles	COLLECT?	3 days before and	If using hormonal supplementation: Discontinue oral and sublingual hormones for 72 hours before and during test. Do not use hormones vaginality, as this may directly contaminate the unine. Other hormones such as topical creams and gels, injectables, and patches can be used on your regular schedule. Do not skip doese of oral birth control pills unless instructed by your healthcare
Men and post-menopausal women Collect any day	Begin the collection between days 19 and 23 of a 28-day cycle, counting the first day of	If cycles are irregular, contact your practitioner for clarity about the timing of	during test	provider. Note: A baseline hormone level may require supplement usage to stop for a longer period of time. Call your health practitioner for more information.
concertary day	your period as day 1.	your collection.	24 hours before and during test	Avoid alcohol, caffeine, tobacco or nicotine-containing products and strenuous exercise. Avoid these foods: avocados, eggplant, tomatoes, bananas, melons, pineapple, grapefruit, plums, fruit juice, nuts, nut butters, wine, cheese, rice, and chocolate.
			Day of Collection	Aveid all supplements and medications until after all samples have been collected (including those that regulate allergy, mood, sleep, pain and inflammation). Consult your practitioner for specific instructions before stopping any medications. Never discontine practicipation medications without first consulting your doctor.
				8 DE DA





# What Are Urinary Hormones and Metabolites?


Unconjugated hormones are sex hormones such as progesterone, testosterone, estradiol, etc.


Conjugated metabolites consist of hormones that the body has modified by adding a water-soluble side arm - usually glucuronide or sulphate molecules - to facilitate their elimination in urine or stool.



























### Can't see the forest for the trees...



• Utilizing urinary hormone testing requires the practitioner to take a step back and look at the pattern of hormone secretion

DCTORS DATA

- Metabolism is complex and will be different for every individual
   Attempting to bring every analyte within range is not necessarily the
- Attempting to bring every analyte within range is not necessarily the goal of this type of testing



### Postmenopausal women who are supplementing


 For postmenopausal women who are supplementing hormones, premenopausal reference ranges are also provided to assist practitioners in assessing treatment

Progesterones		Pre-menopaus Interval	sal Reference							
Progesterone <sup>‡</sup>	(P4)	0.29 – 1.6	Progesterones		Result	Unit	L	WRI	н	Reference Interva
5α-Pregnanediol <sup>‡</sup>	(5A-PD)	50 - 505	Progesterone <sup>‡</sup>	(P4)	0.453	ng/mg Creat/Day			_	0-0.22
Surrregnaneulor	(0410)	50 - 505	5α-Pregnanediol <sup>≇</sup>	(5A-PD)	12.1	ng/mg Creat/Day				21-50
5β-Pregnanediol <sup>‡</sup>	(5B-PD)	375 - 2210	5β-Pregnanediol <sup>‡</sup>	(5B-PD)	121	ng/mg Creat/Day				79 - 280
Allopregnanolone <sup>‡</sup>	(ALLOP)	3.3 – 130	Allopregnanolone <sup>‡</sup>	(ALLOP)	2.55	ng/mg Creat/Day				1.4 - 4.8
04.11		04 50	21-Hydroxyprogesterone <sup>‡</sup>	(21-OHP)	0.639	ng/mg Creat/Day		<b>A</b>		0.3 - 1.4
21-Hydroxyprogesterone <sup>‡</sup>	(21-OHP)	0.4 - 5.6	17-Hydroxyprogesterone <sup>1</sup>	(17-OHP)	0.283	ng/mg Creat/Day				0.17-0.55
5-pregnenetriol <sup>‡</sup>	(5-PT)	70 – 245	5-pregnenetriol <sup>±</sup>	(5-PT)	46.3	ng/mg Creat/Day		A		35 - 120


# Neighborhoods

- Progesterones
- Corticoids
- Androgens
- Estrogens










### Functions of Progesterone

- Opposes effects of estrogen in breast, brain, and endometrial tissue
- Improves estrogen receptor sensitivity
- Essential for maintaining pregnancy
- Promotes sleep patterns (calming effect on the brain/GABA receptors)
- Reduces cholesterol
- Inhibits coronary vasospasm
- Protects against coronary hyper-reactivity
- Diuretic
- Promotes glucose utilization and improves insulin resistance

- Enhances thyroid hormone activity by decreasing thyroid binding globulin
- Increases scalp hair
- Helps to burn fat for energy
- Anti-depressant activity
- Promotes osteoblasts
- Stabilizes and induces secretory endometrial changes
- Promotes cell differentiation
- Promotes normal cell death (apoptosis)

Speroff L. Clinical Gynecologic Endocrinology and Infertility. Philadelphia, PA: Lippincott Williams & Wilkins; 2005. Lee J. Natural Progesterone: The Multiple Roles of a Remarkable Hormone. 4th ed. Sebastopol, CA: BLL Publishing; 1995.



## Progesterone Metabolites

### 5-alpha pregnanediol (5A-PD)

- 5-alpha reduction of progesterone
- Intermediate in the production of allopregnanolone
  - Research indicates 5A-PD may be a partial agonist for GABA receptor-PMID:2550257

### 5-beta pregnanediol (5B-PD)

- Major metabolite of progesterone
- Currently, research indicates this is an inactive metabolite
- Useful as an indirect measurement of progesterone in the body

# Allopregnenolone


- Crosses the BBB binds the GABA receptor
- Research indicates clinical utility for mood disorders, depression, and anxiety
- Paradoxical reaction seen in PMDD-PMID 23978486
- Major metabolite from oral supplementation
- May help evaluate therapeutic oral progesterone metabolism

CHOLESTEROL		P			_		
CVPTAI	PREGNENEDIOL		SI REDUCTANE		ALLOPRE	SNANOLONE	
PREGNENOLONE	FSB282		PROGESTERONE	C7921A2	-> 21-HYDROX	YPROGESTERONE	CATHIBU2
CYPITAL	-> S-PREGNENETRIOL		CYPERAL				
17 OH PREGNENOLONE	F08382	→ 17-HYC	DROXYPROGESTERONE				
· ·							
-							
_							
						<b>D' D</b> 2	



# Corticoids Neighborhood


- Cortisol vs. Cortisone
- Metabolized cortisol
- Cortisol vs. cortisone metabolites



# Cortisol vs. Cortisone


### • Cortisol (active) vs. Cortisone (storage)

- 11B-HSD enzyme: most active in the kidney
  - 11B-HSD1-pulls cortisol out of storage
    - correlated with obesity, metabolic syndrome, inflammation, hypothyroid.
    - 11B-HSD2- protects the mineral corticoid receptor from aldosterone effects
- Because we are looking at the influence of the kidney, we aren't necessarily seeing the HPA axis activity directly, rather the bodies response to stressful stimuli.



# What downregulates 11B-HSD1 activity?

- Reduce inflammation, insulin resistance/insulin, central adiposity
- Physical activity
- Rooibos
- Holy Basil
- Curcumin
- Bitter melon
- EGCG
- Progesterone
- Coffee



# 11B-HSD1 references

- Guilliams T. The Role of Stress and the HPA Axis in Chronic Disease Management. Point Institute; 2015.
- Chapman K, Holmes M, Seckl J. 11β-hydroxysteroid dehydrogenases: intracellular gate-keepers of tissue glucocorticoid action. Physiol Rev. 2013;93(3):1139-1206. doi:10.1152/physrev.00020.2012
- Schloms L, Smith C, Storbeck KH, Marnewick JL, Swart P, Swart AC. Rooibos influences glucocorticoid levels and steroid ratios in vivo and in vitro: a
  natural approach in the management of stress and metabolic disorders?. Mol Nutr Food Res. 2014;58(3):537-549. doi:10.1002/mnfr.201300463
- Jothie Richard E, Illuri R, Bethapudi B, et al. Anti-stress Activity of Ocimum sanctum: Possible Effects on Hypothalamic-Pituitary-Adrenal Axis. Phytother Res. 2016;30(5):805-814. doi:10.1002/ptr.5584
- Hu GX, Lin H, Lian QQ, et al. Curcumin as a potent and selective inhibitor of 11β-hydroxysteroid dehydrogenase 1: improving lipid profiles in high-fattreated rats. PLoS One. 2013;8(3):e49976. doi:10.1371/journal.pone.0049976
- Blum A, Loerz C, Martin HJ, Staab-Weijnitz CA, Maser E. Momordica charantia extract, a herbal remedy for type 2 diabetes, contains a specific 11βhydroxysteroid dehydrogenase type 1 inhibitor. J Steroid Biochem Mol Biol. 2012;128(1-2):51-55. doi:10.1016/j.jsbmb.2011.09.003
- Hintzpeter J, Stapelfeld C, Loerz C, Martin HJ, Maser E. Green tea and one of its constituents, Epigallocatechine-3-gallate, are potent inhibitors of human 11β-hydroxysteroid dehydrogenase type 1. PLoS One. 2014;9(1):e84468. Published 2014 Jan 3. doi:10.1371/journal.pone.0084468
- Atanasov AG, Dzyakanchuk AA, Schweizer RA, Nashev LG, Maurer EM, Odermatt A. Coffee inhibits the reactivation of glucocorticoids by 11betahydroxysteroid dehydrogenase type 1: a glucocorticoid connection in the anti-diabetic action of coffee?. FEBS Lett. 2006;580(17):4081-4085. doi:10.1016/j.febslet.2006.06.046



# Additional corticoid markers


- Metabolized cortisol: (THE + 5αTHF + 5βTHF)
  - (active + potential)
  - measure of what the body has utilized
- Cortisol / Cortisone metabolites: (5αTHF + 5βTHF) vs. THE
  - High metabolized cortisol: increased cortisol clearance due to over production of cortisol, also influenced by; obesity, IR, inflammation, hyperthyroidism
  - Low metabolized cortisol: decreased cortisol clearance which could also be due to hypothyroidism, anorexia, decreased liver function
  - Because free levels of cortisol can be converted to cortisone in the kidney before excretion, looking at the metabolites gives clinicians a better idea of overall metabolic preference within the body






# Androgens Neighborhood

- Evaluate enzyme activity
  - (aromatase and 5 alpha reductase)
- Monitor the activity of the testosterone metabolite 5A-DHT
- 5 alpha vs 5 beta metabolism
  - (Androsterone (5α) / Etiocholanolone (5β))
- DHEA + DHEAs
- Further evaluation of PCOS patients (symptoms, pathways, enzymes)





# Activation of Aromatase (CYP19A1)

- Stress/cortisol: PMID: 21878510 PMID: 23835908 (animal studies, acts in a protective manner for brain health
- Xeno-estrogens/endocrine disruptors/toxins (pesticides, herbicides, benzene, plastic by-products, some pharmaceuticals and cosmetics, petroleum, UV filters) PMID: 2223368
- Poor dietary choices/high glycemic foods PMID: 2223368
- Excess adipose tissue/Obesity/Leptin resistance: PMID: 10349800 PMID: 3226011
- High insulin: PMID: 3322018
- Inflammatory cytokines: IL-6, TNFalpha, prostaglandin PGE(2) PMID: 10405348
- Alcohol/Red wine PMID: 19268535
- Major illness PMID16670151
- Brain injury PMID: 16498364 (animal study)
- Estrogen PMID: 2223368
- Forskolin (found in coleus plant) PMID: 14709151
- Free-fatty acids (increased in metabolic syndrome): PMID: 2223368




# Down regulate Aromatase (CYP19A1)


- Chrysin (honey and propolis): PMID: 32181408
- Zinc PMID 8613886
- Damiana PMID: 18948180
- Resveratrol PMID 16611627
- Plant flavonoids
  - Apigenin PMID 18690828 PMID: 20635310
  - Catechins (EGCG in green tea): PMID 12065214
  - Eriodictyol (Yerba santa) PMID18690828
- Hesperetin (citrus)
- Anti-Mullerian hormone PMID: 2726747
- Glyphosate PMID: 19539684
- Grape seed extract PMID: 16740737
- Nettles: PMID: 17509841
- Ketoconazole PMID: 2004042
- Prolactin

- Metformin PMID: 20300828
- Plant phenols:
  - Chalcones: PMID: 11205867
    - Isoliquiritigenin (licorice) PMID 18690828
- Mangostin PMID 18690828
- Myosime (alkaloid from tobacco and other plants/related to nicotine)
- Nicotine PMID 20188349
- Vitamin E PMID 15084515
- White button mushrooms PMID 17178902
- Aromatase inhibitors (both steroidal and nonsteroidal) i.e. Formestane, Anastrozole
- Things that reduce inflammation:
  - Paclitaxel breast cancer treatment, downregulates TNFreceptors stimulated by aromatase: PMID: 10405348
  - 2-MeOE2- downregulates TNF-receptors stimulated by aromatase: PMID: 10405348
     39



- Enzyme used in the conversion of testosterone to its more potent form; dihydrotestosterone (DHT)
  - DHT has 2-3 times the affinity for the androgen receptor than testosterone
- May lead to increase in androgenic symptoms in females (acne, scalp hair loss, hirsutism, mood instability) and in males, hair loss and prostate issues.






# $\alpha$ vs. $\beta$ : Metabolic Preference of Androgens


Androsterone (5α) /
 Etiocholanolone (5β)

- Alpha (α) metabolism: More potent forms of metabolites
  - May produce more androgenic symptoms
- Beta (β) metabolism: Less potent forms of metabolites
  - May produce less androgenic symptoms




# DHEA + DHEAs


- DHEA is a prohormone with the potential to increase and rogens as well as estrogens
- DHEA is the active form
- DHEAs is the inactive form
- Adding DHEA and DHEAs together may be a better representation of the total level of DHEA than measuring DHEA alone.



DCCTOR'S DATA

43









# Aromatase (CYP19A1)

- Converts androstenedione to estrone and testosterone to estradiol
- Might help explain elevated endogenous estrogens
- Upregulated in peripheral fat in both men and women





# Estrogen Neighborhood

- Aromatase activity-
  - E1 and E2
- Evaluation of metabolites
  - Phase 1 metabolism (-OH)
  - Phase 2 metabolism (-M)
- Evaluation of methylation potential (COMT enzyme)
- Evaluation of potential DNA mutagenesis
- Evaluation of the metabolism of endogenous and exogenous hormones (BHRT)









# CYP1A1

### Upregulate 😳

- DIM/cruciferous veggies PMID: 27261275
- 13C (requires stomach acid to convert to DIM) PMID: 12147290
- Coffee: PMID: 26063478
- Rosemary PMID: 9806165
- Resveratrol PMID: 20716633 PMID: 26167297
- Green and black tea PMID: 24815822
- Fish oil and garlic oil PMID: 12575903
- Hops PMID: 32986415
- Andrographolide, from the Andrographis paniculata plant PMID: 17825862
- Astaxanthin (shrimps and some algae) PMID: 21414371

# CYP1A1

### Downregulate 🗌 😳

- PAHs, PCBs, xenoestrogens, phalates
- BPA PMID: 31374317
- Smoking PMID: 27106177
- Charred meats
- Grapefruit juice
- (bergamottin) PMID: 22705772 • Berries (ellagic acid) can reduce
- overactivity PMID: 26167297 • Green tea extracts PMID: 11064004
- Sulforaphane found in broccoli PMID: 23566952
- St. John's Wort PMID: 16271822
- Lycopene, a red pigment found in tomatoes, carrots, and watermelon PMID: 20400267

- Naringenin and 6',7'-dihydroxybergamottin (from grapefruit juice) PMID: 27444380
- Galangin, found in some plants (Alpinia officinarum, Alpinia galanga, and Helichrysum aureonitens) and propolis PMID: 10188874
- A widely used herbal formulation produced from the extracts of ten common herbs (rosemary, turmeric, ginger, holy basil, green tea, hu zhang, Chinese goldthread, barberry, oregano, and Baikal skullcap) PMID: 22374940




# CYP1B1

- Takes estrogen down the 4 OH pathway
  Pushes E1 to 4-OH E1 and E2 to 4-OH E2
- CYP1B1 is more readily found in tumor tissue compared to normal.
- Can activate cancer promoting compounds
- Accumulating evidence indicates that modulation of CYP1B1 can decrease adipogenesis and tumorigenesis, and prevent obesity, hypertension, atherosclerosis, and cancer. PMID: 28322972
- Given the role of CYP1B1 in pro-carcinogen and estrogen metabolism, polymorphisms in CYP1B1 could result in modifications in its enzyme activity and subsequently lead to hormone-mediated carcinogenesis.
- Also metabolizes: cortisol, aldosterone, FAs, fat-soluble vitamins, melatonin, retinol, plant flavonoids, many environmental toxins.
- Increases fat uptake and can lead to factors of metabolic syndrome

### CYP1B1

### Upregulate 🗌

- Leptin resistance/inflammation/insulin resistance: PMID: 25433128 and DOI: 10.20892/j.issn.2095-3941.2016.0079
- Tetrahydrocannabinol (THC), found in cannabis PMID: 21867498
- UV exposure PMID: 11858729
- Estrogen: PMID 15126349
- PAHs, PCBs: PMID: 28322972
- Diesel exhaust particles (DEP) PMID: 21867498







# CYP1B1

### Downregulate 😳

- Apeacea family: carrots, cumin, anise, celery, carawayPMID: 26381237
- Grapefruit PMID: 16338240
- Resveratrol PMID:32533462
- Quercetin PMID: 22469840 PMID: 16271822
- Apigenin and amentoflavone (St. John's wort) PMID: 16271822
- Ginseng PMID: 11901090
- Lycopene, a red pigment found in tomatoes, carrots, and watermelon PMID: 20400267
- Chrysoeriol, present in rooibos tea and celery PMID: 26167297


- Naringenin (grapefruit juice) PMID: 22935222
- A polyherbal formulation produced from the extracts of ten common herbs (rosemary, turmeric, ginger, holy basil, green tea, hu zhang, Chinese goldthread, barberry, oregano, and Baikal skullcap) PMID:22374940
- Many natural flavonoids and synthetic stilbenes show inhibitory activity toward CYP1B1 expression and function, notably isorhamnetin and 2,4,3',5'tetramethoxystilbene. PMID:28322972
- Genetics:
  - 150 gene polymorphism have been reported



### CYP3A4

### Moves E2 to E3 and E1 to 16OHE1

- It has been estimated that CYP3A4 metabolizes about half of all drugs on the market
- · Mostly found in the liver, but most active in the gut
- The activity of this enzyme varies and is affected by health, environment (smoking, diet, and comedication), hormones, and genetics
- Metabolizes many internal compounds such as cholesterol, fatty acids, prostaglandins, leukotrienes, retinoids and biogenic amines PMID: 25332983.
- Detoxifies bile acids PMID: 25332983
- Partially degrades vitamin D PMID: 22985909



### CYP3A4

### Upregulate

- CYP34A medication inducers
  - Carbamazepine
  - Dexamethasone PMID: 22370628
  - Modafinil
    - Phenobarbital
    - Phenytoin: PMID 10901705
- St.John's wort PMID: 17010103 PMID:1721 4607
- Capsaicin PMID: 22648626
- Valerian PMID: 17214607
- Gingko Bilboa PMID: 17214607

- Vitamin D/UV exposure PMID:22985909
- Being female PMID: 23333322
- Diabetes PMID: 24739263
- Fatty acids PMID: 24739263
- Polycyclic aromatic hydrocarbons (PAH) found in cigarettes PMID: 23845848
- Aflatoxin B1 PMID: 21641981
- Healthy levels of iron (CYP enzymes are heme dependent)




# CYP3A4

### Downregulate

- Polyphenols
  - Flavonoids:
    - Kaempferol (kale, beans, tea, spinach, broccoli) PMID: 25684704
    - Quercetin PMID: 25196644 PMID: 25684704
    - Apigenin PMID: 26180597
    - Chrysin
    - Luteolin
    - Genestein
    - Green tea flavonols EGCG and epicatechin gallate PMID: 26180597
  - Phytoestrogen
    - Coumestrol (soybeans, brussels sprouts, spinach and a variety of legumes)
  - Phenolic acids
    - Caffeic acid PMID: 25196644


- Tannic acid PMID: 26180597
- Gallic acid PMID: 26180597
- Other Polyphenols
  - Licochalcone A from licorice PMID: 26100226
  - Sesamin in sesame seeds (Sesamum indicum) PMID: 22645625
  - Resveratrol PMID: 25341566 PMID: 26180597
  - SulforaphanePMID: 17028159
  - Berberine PMID: 21870106
  - Allyl isothiocyanate (creates the pungent taste of mustard, radish, horseradish, and wasabi) PMID: 25069801
  - Ginsenoside Rd derived from Ginseng PMID: 15133536
  - Gomisin C and gomisin G found in Schisandra Chinensis PMID: 28344076



# Estrogen Metabolism: Phase 2 Detox

### **Oxidation or Methylation?**

- Methylation (COMT)
  - Inactivates catechol estrogens
     Increases solubility/prepare intermediates for renal and biliary excretion
  - methoxy metabolites →are stable and safely leave the
- body Oxidation
  - CYP peroxidase → quinone and semiquinone formation
    - Glutathionization (GST)
    - Glucuronidation (UGT)
    - Sulfonation (SOD2)
    - Acetylation




# Phase 2: Methylation Activity (COMT)

### • Metabolites can demonstrate methylation activity

•The **ratio** of **4-M E1/E2 to 4-OH E1/2** as well as **2-M E1/E2 to 2-OH E1/E2** is evaluated to determine if adequate methylation of catechol estrogens is occurring.

•The higher the ratio of methylation potential, the higher the likelihood of metabolizing more of the catechol estrogens toward the less harmful pathway of methylation and therefore less reactive quinone formation.



Phase 2: Methylation (COMT) Factors

### Upregulate 😳

- Food Support:
  - Cruciferous vegetables
  - Soy foods (ex. genestein)
  - Resveratrol
  - Citrus foods
  - Teas (rooibos, dandelion)
  - Spices (rosemary, curcumin) PMID: 26167297

### **Cofactors and Methyl donors:**

- Methionine
- Vitamin B12
- Vitamin B6
- Betaine
- Folate
- Magnesium

PMID: 26167297

Romilly E. Hodges, Deanna M. Minich, "Modulation of Metabolic Detoxification Pathways Using Foods and Food-Derived Components: A Scientific Review with Clinical Application", Journal of Nutrition and Metabolism, vol. 2015, Article ID 760689, 23 pages, 2015. https://doi.org/10.1155/2015/760689



# Phase 2: Methylation (COMT) Factors

### Downregulate 🗆

- High sucrose diet may inhibit methylation PMID: 26167297
- Leptin resistance: DOI: 10.20892/j.issn.2095-3941.2016.0079
- Estrogen PMID: 10385681 PMID: 30684530
- TNFalpha, present in inflammatory states
- Nuclear factor-kappa B (NF-кB) PMID: 26187567
- Serotonin (competes with the methyl donor S-adenosyl-L-methionine (SAM), thus competing with methylation of COMT) PMID:22500608
- Anything that affects the methionine cycle (homocysteine to cysteine): B6 insufficiency

### Phase 2: Methylation (COMT) Factors

### Downregulate continued . . .

- Bisphenol and PCBs (PMID: 20945454)
- Having too little SAM (s-adenosylmethionine) and too much SAH (sadenosylhomocysteine) from undermethylation
- Entacapone, tolcapone, opicapone and nitecapone
- PPIs, antibiotics
- Genetics polymorphisms PMID: 24593143
  - V158M or rs4680: The A allele results in 3 to 4-fold decrease

# <text><text><text><text><text><text><text><text><text><text>


# Quinone reduction

### NAC

- Reduce estrogen semiquinones back to catechol estrogens
- Primary effect is to react with quinones to form conjugates preventing the formation of estrogen-DNA adducts.

### Resveratrol

- Reduce catechol estrogen semiquinones back to catechol estrogens
- Induce the estrogen-protective enzyme quinone reductase
- Modulates CYP1B1, thereby reducing its activity and thus the formation of 4-OHE<sub>1</sub>(E<sub>2</sub>)
- Shrivastava N, Parikh A, Dewangan RP, et al. Solid Self-Nano Emulsifying Nanoplatform Loaded with Tamoxifen and Resveratrol for Treatment of Breast Cancer. *Pharmaceutics*. 2022;14(7):1486. Published 2022 Jul 18. doi:10.3390/pharmaceutics14071486



# 8-hydroxy-2-deoxyguanosine (80H2dG)

### · Not an estrogen metabolite, but pairs nicely with metabolite information

- Marker for DNA damage due to oxidative stress in general
- May be helpful in identifying or confirming DNA damage from harmful metabolites (quinones and DNA adducts)
- In nuclear and mitochondrial DNA, 8-OHdG is one of the predominant forms of free radical-induced oxidative lesions and widely used as a biomarker for oxidative stress and carcinogenesis.

8-hydroxy-2'-deoxyguanosine <sup>4</sup> (8-OHdG) 9.08 ng/mg Creat/Day 0 – 7.5	Oxidative Stress Metabolite	Result	Unit	L	WRI	1	н	Refere	nce Interval
	8-hydroxy-2'-deoxyguanosine <sup>‡</sup> (8-OHdG)	9.08	ng/mg Creat/Day					0-7.5	

### Oxidative stress and disease states Chronic stress Cystic fibrosis Cortisol elevation Psoriasis Inflammation Chronic hepatitis Insomnia in postmenopausal women Gastritis • Degenerative diseases (rheumatoid arthritis, Irritable boweldisease Parkinson's disease, Huntington's disease, Pancreatitis Alzheimer's disease) Cancer Chronic fatigue syndrome Acute viral infection Major depression • Copper implants Hypertension · Toxic exposures (tobacco smoke, Cardiovascular disease methamphetamines, asbestos, heavy metals, Diabetes type II polycyclic hydrocarbons) DOCTOR'S DATA SCIENCE + INSIGHT

### Treatment for elevated 8-OHdG

### Address the cause of oxidative stress

- Glutathione
- NAC
- Green tea
- CoQ10
- Onion
- Garlic
- · Alpha lipoic acid
- Vitamins C and E

- Melatonin
- Folate
- Berberine
- EPA/DHA
- Fermented papaya powder
- Increase fruit and veggies
- Yoga



### 8-OHdG:treatment references


- Semenova NV, Madaeva IM, Brichagina AS, Kolesnikov SI, Kolesnikova LI. 8-Hydroxy-2'-Deoxyguanosine as an Oxidative Stress Marker in Insomnia. Bull Exp Biol Med. 2021;171(3):384-387. doi:10.1007/s10517-021-05233-0
- Sidorova Y, Domanskyi A. Detecting Oxidative Stress Biomarkers in Neurodegenerative Disease Models and Patients. Methods Protoc. 2020;3(4):66. Published 2020 Sep 24. doi:10.3390/mps3040066
   BarbagalloM, Marotta F, Dominguez LJ. Oxidative stress in patients with Alzheimer's disease: effect of extracts of fermented papaya pow der. Mediators Inflamm. 2015;2015;624801. doi:10.1155/2015/624801
- Ceylan D, YilmazS, Tuna G, et al. Alterations in levels of 8-Oxo-2<sup>+</sup>-deoxyguanosine and 8-Oxoguanine DNA glycosylase 1 during a current episode and after remission in unipolar and bipolar depression. *Psychoneuroendocrinology*, 2020;114:104600. doi:10.1016/j.psyneure.2020.104600
- Di Minno A, Turnu L, PorroB, et al. & Hydroxy-2-Deoxyguanosine Levels and Cardiovascular Disease: A Systematic Review and Meta-Analysis of the Literature. Antioxid Redox Signal. 2016;24(10):548-555. doi:10.1089/ars.2015.6508
- Ye X, Jiang R, Zhang Q, et al. Increased 8-hydroxy-2'-deoxyguanosine in leukocyte DNA from patients with type 2 diabetes and microangiopathy. JInt Med Res. 2016;44(3):472-482.doi:10.1177/0300060515521530
- Nakai K, Yoneda K, Maeda R, et al. Urinary biomarker of oxidative stress in patients with psoriasis vulgaris and atopic dermatitis. J Eur Acad Dermatol Venereol. 2009;23(12):1405-1408. doi:10.1111/j.1468-3083.2009.03327.x
- Ock CY, Kim EH, Choi DJ, Lee HJ, Hahm KB, Chung MH. 8-Hydroxydeoxyguanosine: not mere biomarker for oxidative stress, but remedy for oxidative stress-implicated gastrointestinal diseases. World J Gastroenterol. 2012;18(4):302-308. doi:10.3748/wjg.v18.i4.302
- Jelic MD, Mandic AD, Maricic SM, Srdjenovic BU. Oxidative stress and its role in cancer. J Cancer Res Ther. 2021;17(1):22-28. doi:10.4103/jcrt.JCRT\_862\_16
- Toyokuni S, Sagripanti JL. Increased 8-hydroxydeoxyguanosine in kidney and liver of rats continuously exposed to copper. Toxicol Appl Pharmacol. 1994;126(1):91-97. doi:10.1006/taap.1994.1094
- Huang MC, Lai YC, Lin SK, Chen GH. Increased blood 8-hydroxy-2-deoxyguanosine levels in methamphetamine users during early abstinence. Am J Drug Alcohol Abuse. 2018;44(3):385-402. doi:10.1080/00952990.2017.1344683
- Valavanidis A, Vlachogianni T, Flotakis C. 8-hydroxy-2'-deoxyguanosine (8-OHdG): A critical biomarker of oxidative stress and carcinogenesis. J Environ Sci Health C Environ Ecotoxicol Rev. 2009;27(2):120-139. doi:10.1080/1059050090288568
- Arab H, Mahjoub S, Hajian-Tilaki K, Moghadasi M. The effect of green tea consumption on oxidative stress markers and cognitive function in patients with Alzheimer's disease: A prospective intervention study. Cospion J Intern Med. 2016;7(3):188-194.
- Biglan KM, Dorsey ER, Evans RV, et al. Plasma &-hydroxy-2'-deoxyguanosine Levels in Huntington Disease and Healthy Controls Treated with Coenzyme Q10. J Huntingtons Dis. 2012;1(1):65-69. doi:10.3233/JHD-2012-120007
- Boyle SP, Dobson VL, Duthie SJ, Kyle JA, Collins AR. Absorption and DNA protective effects of flavonoid glycosides from an onion meal. Eur J Nutr. 2000;39(5):213-223. doi:10.1007/s003940070014
- Zalejska-Fiolka J, Wielkoszyński T, Rokicki W Jr, et al. The Influence of α-Lipoic Acid and Garlic Administration on Biomarkers of Oxidative Stress and Inflammation in Rabbits Exposed to Oxidized Nutrition Oils. Biomed Res Int. 2015;2015;827879. doi:10.1155/2015/827879
- Guest J, Bilgin A, Hokin B, Mori TA, Croft KD, Grant R. Novel relationships between B12, folate and markers of inflammation, oxidative stress and NAD(H) levels, systemically and in the CNSof a healthy human cohort. Nutr Neurosci. 2015;18(8):355-364. doi:10.1179/1476830515Y.0000000041
- Dai P, Wang J, Lin L, Zhang Y, Wang Z. Renoprotective effects of berberine as adjuvant therapy for hypertensive patients with type 2 diabetes mellitus: Evaluation via biochemical markers and color doppler ultrasonography. Exp Ther Med 2015, Sep. 10(3):869-76.
- Ghorbanihaghjo A, Safa J, Alizadeh S, et al. Protective effect of fishoil supplementation on DNA damage induced by cigarette smoking. J Health Popul Nutr. 2013;31(3):343-349. doi:10.3329/jhpn.v3113.16826
- Barbagallo M, Marotta F, Dominguez LJ. Oxidative stress in patients with Alzheimer's disease: effect of extracts of fermented papaya pow der. Mediators Inflamm. 2015;2015:624801.doi:10.1155/2015/624801
- Franke AA, Cooney RV, Henning SM, Custer LJ. Bioavailability and antioxidant effects of orange juice components in humans. JAgric Food Chem. 2005;53(13):5170-5178.doi:10.1021/jf050054y
- Nirwan M, Halder K, Saha M, Pathak A, Balakrishnan R, Ganju L. Improvement in resilience and stress-related blood markers following ten months yoga practice in Antarctica. J Complement Integr Med. 2020;18(1):201-207. Published 2020 Jun 19. doi:10.1515/jcim-2019-0240

#### What are the advantages of HuMap<sup>™</sup>?

#### **Additional Analytes**

- Progesterone itself
- DHEA
- Allopregnanolone
- Androstenedione
- (2-OH E2 and 4-OH E2) included in our ratios
- (2-M E2 and 4-M E2) included in our ratios


#### **Liquid Collection**




• With liquid urine, samples can be shipped after being frozen for 4-6 hours, can be processed faster, and concentrated further to enhance the detection of low-level analytes. Steroid hormones and metabolites are also quite stable in liquid urine if the correct preservative is used.

metabolites or creatinine for some patient samples.





#### Elimination or Recirculation




Baker, J. M., Al-Nakkash, L., & Herbst-Kralovetz, M. M. (2017). Estrogen-gut microbiome axis: Physiological and clinical implications. *Maturitas* 103, 45-53. doi:10.1016/j.maturitas.2017.06.025

- Conjugated estrogens pass into the gut for elimination
- Beta-glucuronidase can deconjugate estrogens
- Unconjugated estrogens reabsorbed into circulation
- High levels of beta-glucuronidase could mean more recirculating of estrogens
- Beta-glucuronidase increases with: low fruits and veggies, low fiber, high sugar, processed foods, SAD diet, alcohol, toxicants, antibiotics

DCTORS DATA

50



#### Consider treating phase 3 (elimination) first

- Hydration PMID: 20646222
- Fiber PMID: 26026145
- Herbs to prevent/treat constipation (aloe, ginger, Avipattikar, Triphala) PMID: 30680163 PMID: 28696777
- Magnesium
- Movement
- Calcium D-Glucarate PMID: 2346674

Genuis SJ. Elimination of persistent toxicants from the human body. Hum Exp Toxicol. 2011 Jan;30(1):3-18. doi: 10.1177/0960327110368417. Epub 2010 Apr 16. PMID: 20400489.

#### Urinary Hormone and Metabolites Summary


- HuMap<sup>™</sup> gives a comprehensive assessment of unconjugated hormones and their metabolites
- Allows the practitioner to examine what the body has utilized
- Can give greater insight into
  - Alpha vs. Beta metabolism (progesterones, androgens, and corticoids)
  - Phase I and Phase II detoxification (COMT /methylation activity)
  - With additional estradiol metabolites for 2-OH and 4-OH to give a greater picture of metabolism
  - Potential risk assessment for breast health and other disease processes
  - Utilization of hormones both endogenous and exogenous
  - Can assist with complex cases where BHRT and salivary testing appear to be normal, but the patient
    may still have symptoms.
- Remember, elimination is key to any treatment case
  - Ensure proper elimination first before addressing phase 1 or phase 2 metabolism



## Why test neurotransmitters and urinary metabolites together?

It could be especially important to combine HuMap<sup>™</sup> and NT testing in cases involving mental health symptoms that came on during a time of hormonal transition:


- For example, an aging male with depression might have low testosterone.
- Mood symptoms that happen at certain times of the menstrual cycle can involve NT imbalance and also estrogen dominance.





DCCTORS DATA


#### COMT

• An obvious connection between HuMap<sup>™</sup> and Neurotransmitter testing is COMT activity, which is important with estrogen and catecholamine metabolism.









#### **COMT** Support

#### **COMT** sluggish

- SAMe: 100-500 mg
- Mg: 150-500 mg
- MTHF: 400-5000 mcg
- Methylcobalamin: 1000-5000 mcg


Davis 2020.


- Foods/herbs: cruciferous vegetables, soy foods (ex. genestein), resveratrol, citrus, rooibos, dandelion, rosemary, curcumin
- Limit caffeine, alcohol, tobacco smoke

#### **COMT** over active

- Support MAOenzyme...why?
  - Vitamin B2: 50 mg (riboflavin 5 phosphate)

  - Vitamin B3: 100 mg
    - Iron: 25-50 mg






#### Other factors that affect enzyme activity

- Inflammatory cytokines can influence the activity of both MAO and COMT.
- Insulin resistance may increase MAO levels
- Excess estradiol slows COMT and MAO activity, potentially resulting in higher levels of neurotransmitters
- Low estradiol increases MAO activity, lowering serotonin
- Progesterone downregulates COMT gene expression

Salama, Salama A et al. "Progesterone regulates catechol-O-methyl transferase gene expression in breast cancer cells: distinct effect of progesterone receptor isoforms." The Journal of steroid biochemistry and molecular biology vol. 107,3-5 (2007): 253-61. doi:10.1016/j.jsbmb.2007.03.049 Kleinridders A, Cai W, Cappellucci L, et al. Insulin resistance in brain alters dopamine turnover and causes behavioral disorders. Proc Natl Acad Sci U S A. 2015;112(11):3463-3468. doi:10.1073/pnas.1500877112



DOCTORS DATA

#### Coordinated fight or flight response

- Norepinephrine and epinephrine are released from the brain as well as storage vesicles in the adrenal medulla in response to stress: fright, exercise, cold, low blood glucose.
- They increase degradation of glycogen and triacylglycerol, as well as increase blood pressure and cardiac output.
- Cortisol is then released to continue the stress response and continues to be released until the threat is over.



Champe P, Harvey R, Ferrier D. Biochemistry. 3rd ed. Philadelphia: Lippincott Williams & Wilkins; 2005.

#### Stress and the Brain

- Stress literally shrinks the brain and degrades the BBB
  - Stressors: Environmental stress, but also smoking, food intolerances, blood sugar imbalances, anemia, bacterial gut infections, gut parasites, autoimmune dz, joint pain and inflammation, poor digestion, etc.
- Excess fat is pro-inflammatory and a chronic stressor for body and brain.
- Studies show that high cortisol in response to high stress damages the hippocampus, which regulates our circadian rhythm.
  - Eventually the adrenals become less responsive to ACTH.


Esposito P, et al. Acute stress increases permeability of the blood-brain-barrier through activation of brain mast cells. Brain Res. 2001 Jan5;888(1):117-127. Tavanti M et al. Evidence of diffuse damage in frontal and occipital cortex in the brain of patients with post-traumatic stress disorder. NeurolSci. 2012 Feb;33(1):59-68.

#### Stress response and serotonin


Stress response


• The response of the CNS to stress is predominantly through activation of the locus coeruleus/sympathetic nervous system, with the consequent release of **norepinephrine and epinephrine**, and the limbic/HPA system, with the consecutive release of **serotonin**, CRH, ACTH, and **cortisol**.


Tafet GE, Idoyaga-Vargas VP, Abulafi a DP, et al. Correlation between cortisol level and serotonin uptake in patients with chronic stress and depression. Cogn Affect Behav Neurosci. 2001 Dec; 1(4): 388-93.





# <section-header><list-item><list-item><list-item><list-item><list-item><list-item><list-item><table-container>



























Ratios and Calculation	s	Result	Unit	L	WRI	Н	Reference Interval
2-OH-E1 % <sup>‡</sup>	(2-OH-E1 %)	60.2	%				50 – 85
4-OH-E1 % <sup>‡</sup>	(4-OH-E1 %)	3.23	%				2-7
16-OH-E1 % <sup>‡</sup>	(16-OH-E1 %)	36.6	%				15 – 50
2-M-E1:2-OH-E1 <sup>‡</sup>	(COMT/Methylation activity)	0.395					0.1-0.36
2-M-E2:2-OH-E2 <sup>‡</sup>	(COMT/Methylation activity)	0.348					0.07 - 0.37
4-M-E1:4-OH-E1 <sup>‡</sup>	(COMT/Methylation activity)	0.195					0.09-0.54
4-M-E2:4-OH-E2 <sup>‡</sup>	(COMT/Methylation activity)	0.095					0.04 - 0.54
2-OH-E1:16-OH-E1 *		1.64					1.6 – 5.1
4-OH-E1:2-OH-E1 <sup>‡</sup>		0.054					0.02-0.07
Oxidative Stress Metab	polite	Result	Unit			. — — —	Reference Interval
8-hydroxy-2'-deoxyguan	osine <sup>‡</sup> (8-OHdG)	9.20	ng/mg Creat/Day				0-7.5




#### Thyroid

- Thyroid issues can affect both NT secretion and urinary metabolite secretion.
- HuMap<sup>™</sup> Hypothyroid could result in higher cortisol, low cortisol metabolites.
  - Lower metabolites in general, across the report (slower metabolism)
  - Hyperthyroid could result in high cortisol metabolites, elevated cortisone.

#### • NT testing

- Low tyrosine could impact thyroid.
- Hypothyroid = low dopamine, low serotonin, low GABA, and could increase toxicity of elevated glutamate.



5A-THF+5B-THF/THE <sup>‡</sup>	(Cortisol/Cortisone Metabolites)	0.832	<b>A</b>		0.7 – 1.4
Cortisol/Cortisone <sup>‡</sup>	(11B HSD activity)	0.638			0.24 - 0
	CORTISOL HSOTHEC/HSOTHE CORTISONE	-	B-THF THE		
				107	DOCTORS DATA


Cortisol/Cortisone <sup>±</sup>	(11B HSD activity)	0.184	A		0.24-0.45
5A-THF/5B-THF ratio*	(alpha vs beta metabolism)	0.902		_	0.4-1.4
	CORTISOL HSDTTE2 / HSDTTE1 CORTISONE Possibly ass	ociate	SG REDUCTASE / HSD3A $\rightarrow$ 5A-THF SG REDUCTASE / HSD3A $\rightarrow$ 5B-THF SG REDUCTASE / HSD3A $\rightarrow$ THE d with hyperthyroidism		
				108	DCTIONE DIAL


Analyte	Result	Unit per Creatinine	L	WRI	н	Reference Interval	
Phenethylamine (PEA)	20	nmol/g				42-160	
Tyrosine	48	µmol/g				70–180	
Tyramine	3.4	µmol/g	_			2.8-8.5	
Dopamine	258	hð\ð				175 – 500	
3,4-Dihydroxyphenylacetic acid (DOPAC)	1520	hð\ð				540 - 1850	
3-Methoxytyramine (3-MT)	151	nmol/g				122-278	
Norepinephrine	17.6	hð\ð				29-69	
Normetanephrine	195	hð\ð				112-400	
Epinephrine	1.9	hð\ð				2.1-14.5	
Metanephrine	84	hð\ð				60-158	
Norepinephrine / Epinephrine ratio	9.3					< 13	
Tryptamine	0.4	µmol/g				0.65 - 1.6	
Serotonin	84.4	hð\ð				79–235	
5-Hydroxyindoleacetic acid (5-HIAA)	5913	hð\ð				2500-9000	
Glutamate	23	nmol/g		<b>A</b>		18.0-70.0	
Gamma-aminobutyrate (GABA)	5	nmol/g				2.6-8.0	
Glycine	966	nmol/g				700-2500	
Histamine	12	hð\ð				14 – 51	
Taurine	428	µmol/g	Δ			420-1400	
Creatinine	73.1	mg/dL				25-180	

#### Sex Hormones affect NT secretion and vice versa

- Estradiol enhances serotonin and dopamine levels in males and females
- Progesterone (allopregnanolone) stimulates GABA receptors in males and females
- Dopamine enhances testosterone secretion in males (bidirectional)
- NE and insulin are increased in metabolic syndrome. Insulin drives the metabolic changes and NE drives the hypertension and impaired glucose metabolism.










### Estradiol increases serotonin synthesis and decreases its breakdown


- Higher levels of estradiol inhibit MAO leading to more serotonin available in the synapse.
- Estradiol supplementation increases tryptophan hydroxylase
- Increased serotonin activity improves cognitive performance and mood.

Amin Z, Canli T, Epperson CN. Effect of estrogen-serotonin interactions on mood and cognition. *Behav Cogn Neurosci Rev.* 2005; 4: 43-58. Epperson CN, Amin Z, Ruparel K, et al. Interactive effects of estrogen and serotonin on brain activation during working memory and affective processing in menopausal women. *Psychoneuroendocrinology.* 2012; 37(3): 372-382 Klaiber EL, Broverman DM, Vogel W, et al. Effects of estrogen therapy on plasma MAO activity and EEG driving responses of depressed women. Am J Psychiatry. 1972 Jun; 128(12): 1492-8. Paredes S, Cantillo S, Candido KD, Knezevic MN. An Association of Serotonin with Pain Disorders and Its Modulation by Estrogens. *Int J Mol Sci.* 2019;20(22):5729. Published 2019 Nov 15.









# Dopamine and testosterone in malelibido Relationship is bidirectional; dopamine can influence testosterone and testosterone can influence dopamine. In males, an important area of thebrain for sexual function is the medial preoptic area (MPA). One study found that microinjecting dopamine agonists in the MPA of rats resulted in an increase of sexual activity Another study observed that castrated male rats showed no interest in

• Another study observed that castrated male rats showed no interest in fornicating and no dopamine release in the medial preoptic area. After testosterone injections, the castrated rats were able to engage in sexual intercourse and showed an increase in dopamine release in the MPA.

Dominguez JM, Hull EM. Dopamine, the medial preoptic area, and male sexual behavior. *Physiol Behav* 2005;86: 356-368 Putnam SK, Du J, Sato S, Hull EM. Testosterone Restoration of Copulatory Behavior Correlates with Medial Preoptic Dopamine Release in Castrated Male Rats. *Hormones and Behavior* 2001;39(3):216–224.

#### Sex Hormone and NTs Summary

- Estradiol enhances serotonin and dopamine levels in males and females
- Progesterone (allopregnanolone) stimulates GABA receptors in males and females
- Dopamine enhances testosterone secretion in males (bidirectional)



#### Independent Influence of Insulin and Catecholamines on Metabolic Syndrome

- Metabolic syndrome is a state of sympathetic nerve hyperactivity.
  - Associated with elevations in NE and epi, but mostly NE.
  - Insulin and NE levels increased in Met Synd.
- This study suggests that insulin and NE cooperate independently to the development of metabolic syndrome.

De Pergola G, Giorgino F, Benigno R, et al. Independent influence of insulin, catecholamines, and thyroid hormones on metabolic syndrome. Obesity (SilverSpring). 2008 Nov; 16(11): 2405-11.

#### Independent Influence of Insulin and Catecholamines on Metabolic Syndrome

- Insulin and catecholamines likely cooperate to the development and clinical features of met syndrome.
- Insulin seems to be the determinant for all the negative metabolic aspects, whereas catecholamines seem to be important for the development of hypertension and impaired glucose metabolism.



De Pergola G, Giorgino F, Benigno R, et al. Independent influence of insulin, catecholamines, and thyroid hormones on metabolic syndrome. Obesity (Silver Spring). 2008 Nov; 16(11): 2405-11.

#### Symptomology

- And if we just think about symptomology, working up hot flashes would benefit from both tests, as they are associated with estrogen, progesterone and serotonin imbalance.
- Anxiety, depression, and cognitive symptoms are associated with estrogen dominance, low progesterone, low estradiol, as well as various NT imbalances.
- Not all symptoms result in predictable patterns, so if there is the possibility of there being a neuroendocrine component to a symptom or pathology, test both!

Patient: Sample Patient         Sample Collection Date/Time           Mit 200909.0000.1         Age:01 DOB:0101/1960         Sample Collection Date/Time           Bedtime         12/30/2022 22.30	Order: SAMPLE REPORT Patient: Sample Patient Sample Collection Date/Time
kert #: 199999     kert #: 19999     kert #: 19999     kert #: 1999     kert #: 1999     kert #: 1991	Distance
The Barrymph regression the inflationality of the assession activations 12 (344) 4 (344), 14 (3	Analytica Beauty Linit per Creatinine L WRI H Reference Interval
Properture	nolig 42-180
	aulia 70.180
B) Propriedbi	1919 175-500
	19/9 A 1850
Property Dista-CHAS	nolig 122-278
	ig/g 29-69
Continential	igig 112-400
0	19/g 2.1-14.5
Hetabolised Certisel     cartitat to confidence. Inhibition of this engine     may lead to the amount of control heing	Metanephrine 84 µg/g 60-158
greater films configure, while increased entypes activity can base to higher levels of	Norepinephrine / Epinephrine ratio 9.3 <
	Tryptamine 0.4 μmol/g 0.65-1.6
	Serotonin 84.4 μg/g 79-235
C Allopageaulore Central Contrare	5-Hydroxyindolescetic acid (5-HIAA) 5913 µg/g 2500 - 9000
TireS21 Weng Tri Fait Weng Browting Browting	Giutamate 23 nmol/g 18.0-70.0
KEY BELATIONSHIPS     KEY BELATIONSHIPS     To spuids to the non-invalidation     EA REDUICTISSE ACTIVITY	Gamma-aminobutyrate (GABA) 5 nmolig 2.6-8.0
preference by key stappen, indicated by the	Giycine 966 nmolig 700-2500
Piofabulhas in the back pathway are more provided in the pione of the back pathway are more provided in the back provided in the pione of the back provided in the pione of th	Histamine 12 µg/g 14-51
Cash be responsible for androgenic services and the second services and the services and the second se	Taurine 428 µmolig 420 – 1400
Accurations is an ensyme found in the groutest and an ensyme found in the groutest and an ensyme found in the groutest and the ensyme ensymm ensymmetry ensym	Creatinine 73.1 mg/dL 25 - 180
Constanting of the second and the se	Neurotransmitter Comments: Urlang neurotransmitter levels provide an overall assessment of the body's ability to make and break down neurotransmitters and are recommentative of which body levels. Neurotransmitters are secreted all through the body, in neurons of both the control and peripheral nervous
P 1-Presented	systems. The express, collectors and programmers in neurotensmitter metabolism in prevent are the same in the pumpting valid in the central memory system. Therefore, that laterotary is unary monotonemers there is a same set of the pumptant clinical information, and may be associated with many symptoms includery cognitive and mode concerns, diminished drive, futigo and steep diffusions, carvings, addictions and a set of the metabolism includers of the dimensional set of the d
4-HYDROLEVED HOME WE AMERICAN HIS HITME	is the precursor arrive call to PEA, and vitamin B6 is a required co-factor in the conversion to this primary trace arrive.