

Introduction to Epigenetic Methylation Biomarkers

About TruDiagnostic: Epigenetic Methylation Specialists



-  Over 35 institutional research board approved clinical trials
-  20,000 sq ft CLIA Lab 3 Custom Aging and Ingestions Kits
-  Started in July of 2020. Have tested over 10,000 patients in that time with multiple operators
-  Collaborative Research with the top epigenetic research institutions

TruDiagnostic™

What is Epigenetics?




What is Epigenetics?

Every cell in your body has the exact same DNA sequence.

So how do your cells act so differently?


It is due to what DNA expression is turned on and turn off!




What is Epigenetics?

Epigenetics is


DNA Methylation



Histone Modification



Non-coding RNA



"the branch of biology which studies the causal interactions between genes and their products which bring the phenotype into being."
- Conrad Waddington 1942s

TruDiagnostic


What is DNA Methylation?


DNA methylation is a biological process by which methyl groups are added to the DNA molecule. Methylation can change the activity of a DNA segment without changing the sequence.

When located in a gene promoter, DNA methylation typically acts to repress gene transcription

There are approximately 29 million locations in the genome where we can methylate and it is different in every cell.

Typical mammalian DNA methylation landscape





What do we measure? How much is measured?

Definition of Beta-value and M-value
 The Beta-value is the ratio of the methylated probe intensity and the overall intensity (sum of methylated and unmethylated probe intensities). Following the notation used by Illumina methylation assay (1), Beta-value for an I^{th} interrogated CpG site is defined as:

$$\text{Beta}_i = \frac{\text{max}(Y_{i,\text{methyl}}, 0)}{\text{max}(Y_{i,\text{methyl}}, 0) + \text{max}(Y_{i,\text{unmethyl}}, 0) + \alpha} \quad (1)$$


where $Y_{i,\text{methyl}}$ and $Y_{i,\text{unmethyl}}$ are the intensities measured by the I^{th} methylated and unmethylated probes, respectively. To avoid negative values after background adjustment, any negative values will be reset to 0. Illu-



TruDiagnostic™

Why should we measure? The Famous Agouti Mouse Experiment


These two mice are Genetically Identical and the Same Age



While pregnant, both of their mothers were fed Bisphenol A (BPA) but DIFFERENT DIETS:


The mother of this mouse received a normal mouse diet.
 The mother of this mouse received a diet supplemented with choline, folate, zinc, betaine and vitamin B12.

How Methylation Can Influence Obesity: The Agouti Gene



TruDiagnostic™

Methylation as a Biomarker and Treatment Target




Biomarkers have been defined as: "indicators of biological and pathogenic processes, or pharmacologic responses to a therapeutic intervention that defines what is normal while predicting or detecting what is abnormal."

TruDiagnostic™

Methylation as a Biomarker

5 to 10 years from now, the health system that doesn't use this data to improve their medical delivery is going to be deemed archaic.



Peripheral blood DNA methylation-based machine learning models for prediction of knee osteoarthritis progression: bioinformatics and data from the Osteoarthritis Initiative and Johnston County Osteoarthritis Project

Journal of Bone and Joint Surgery | Volume 94-B | July 2012 | 1313-1319
DOI: 10.1302/0301-2203.94B07.2011.0301161

Clinical and Translational Medicine

Comprehensive methylation sequencing reveals prognostic epigenetic biomarkers for prostate cancer mortality


Journal of Clinical Oncology | Volume 30, Number 4 | February 1, 2012 | Pages 451-458


Methylation risk scores are associated with a collection of phenotypes within electronic health record systems

Journal of the American Medical Association | Volume 305, Number 12 | March 21, 2011 | Pages 1338-1344


TruDiagnostic

Measurement means nothing without interpretation of the data





Linear Model
Logistic Model
Neural Network Model



Artificial Intelligence


- Ability of a machine to make intelligent human decisions
- Application of AI that allows a system to automatically learn and improve from experience
- Application of Machine Learning that uses complex algorithms and neural networks to make a model

TruDiagnostic

What can Epigenetic Measurements Tell Us?

Tru's epigenetic methylation methodology can be trained to predict all other multi-omic values

The Multiome



Methylation is different in every cell and each cell has over 26 million CpG methylation locations.

Due to its large coverage and its place in the hierarchy of the multi-ome, it is able to be trained more accurately than most other multi-omic biomarkers to predict many different outcomes with a single test.

Beyond disease and phenotypes, methylation markers have been shown to predict prior exposures/lifestyle, future outcomes, and current behaviors

Organization (genetic, cell, organ, organism)

Age/Environment

Disruption of Multi-Omic Data

TruDiagnostic

Possibly, ONE test ... to replace MANY

One test to predict many chronic diseases, exposures and other omic biomarkers

TruDiagnostic becomes the best aging and disease predictor available on the market

Methylation trained to predict **health phenotypes** (diagnostics, prognosis, status):

- Efficacy and side effects response to medication
- Death prediction
- Biological aging rate
- Schizophrenia diagnosis
- Cancer (multiple types at as early as stage zero)
- Viral infections
- Immunodeficiencies

Methylation trained to predict **other omic biomarkers**:

- DNA epigenetic acetylation markers
- Senescence
- Adrenomedullin
- Beta-2-microglobulin
- CDS5
- Cystatin C
- EGF ECM protein 1
- Legtin
- Myoglobin
- Plasminogen activator inhibitor
- Tissue inhibitor metalloproteinase
- Immune cell subsets
- Senescent T cells
- SNP polymorphisms

Methylation trained to predict **exposures**:

- Alkaloids B1
- Air pollution
- Arsenic
- Bisphenol A
- Cadmium
- Chromium
- Lead
- Mercury
- Polycyclic aromatic hydrocarbons
- Persistent organic pollutants
- Smoking status (active, former, current)
- Smoking pack years
- Current alcohol intake
- Radiation exposure
- Nutrient status

TruDiagnostic™

Arthritis

Cognitive Decline

Immune Decline

Type II Diabetes

Autoimmune Disorders

Cancer

Kidney Disease

Cardiovascular Disease

?

What is the #1 risk factor for all of these conditions?

TruDiagnostic™

TOP 3 CAUSES OF DEATH FOR 2020 IN THE U.S.

Heart Disease


Cancer (all)

COVID-19

?

What is the #1 risk factor for all of these conditions?

TruDiagnostic™



Aging As a Disease

Aging is the **#1 risk factor** for almost ALL chronic diseases.

However, delaying the aging rate by 7 years would **cut the incidence of disease in half!**

If everyone were to increase life expectancy by 1 year, we would save **\$256.7 billion** and **\$8.8 trillion** for +10 years life expectancy.

TruDiagnostic™

Since aging is the greatest risk factor for every disease, how can we appropriately measure it? Are there better methods than chronological age?


- Chronological age (CA) is a commonly used indicator for aging. However, life expectancy shows considerable variation among individuals with equal or similar CAs due to diversity in genotypes and in living habits and environments.
- A 50-year-old individual may have 60-year-old body functions, and many people look older or younger compared to others at the same CA (even in twins). Therefore, CA is not an optimal indicator for the aging progress.

The History Of Biological Age Measurements

During the past decades, extensive effort has been made to identify such aging biomarkers that, according to the stage-setting definition (Baker and Sprott, 1988), are "biological parameters of an organism that either alone or in some multivariate composite will, in the absence of disease, better predict functional capability at some late age, than will chronological age". Later on, the American Federation for Aging Research (AFAR) formulated the criteria for aging biomarkers as follows:

1. It must predict the rate of aging. In other words, it would tell exactly where a person is in their total life span. It must be a better predictor of life span than chronological age.
2. It must monitor a basic process that underlies the aging process, not the effects of disease.
3. It must be able to be tested repeatedly without harming the person. For example, a blood test or an imaging technique.
4. It must be something that works in humans and in laboratory animals, such as mice. This is so that it can be tested in lab animals before being validated in humans.

However, until 2013, no such marker or marker combination had emerged.



How does this compare to telomeres?

"Briefly, telomere length is extensively validated but has low predictive power."

EBioMedicine

Journal Pre-proof

Biological Age Predictors
David A. Huffman, Nancy L. Pedersen, Sara Klipp

ARTICLE INFO

ABSTRACT

We need a way to track age... Thankfully we finally have one. Epigenetics Clocks!

5 to 10 years from now, the health system that doesn't use this data to improve their medical delivery is going to be deemed archaic.

- Atul Butta, President of TruDiagnosis, Researcher in Genetic Testing



TruDiagnostic

How should someone use Epigenetic Age Testing?

- Epigenetic Age = Biggest Risk Factor For Chronic Disease and Death
- Epigenetic Age = Major correlate of quality of life metrics like IQ, Balance, Muscle Mass, and Facial Aging
- Accelerated Aging = Shorter Lifespan and Healthspan
- Decelerated aging = Longer Lifespan and Healthspan
- Aging starts as soon as we are born. So aging optimization early can be a great preventative tool.

Why wouldn't everyone try and measure and treat their major risk factor for all disease? (especially if they are currently "healthy")

The Predictive Power of Epigenetic Aging: Cancer

- ◆ Cancer is an age-related condition linked to epigenetic age acceleration
- ◇ "For each one year increase in the differences between chronological and epigenetic age (the Δage), there was a 6% increased risk of developing cancer within three years and a 17% increased risk of dying of cancer in the next five years" (Jylhävä et al 2017)
- ◆ "Some cancers, tumours with fewer somatic mutations, indicative of more stable genomes, have greater epigenetic age acceleration" (Horvath & Raj 2018)

TruDiagnostic's Focus on Aging

Patient Reports

TruAge
Built for contextual benefit to the patient and for epidemiological review on epigenetic aging

Pace Of Aging
One of the most accurate, precise and most predictive of age-related outcomes

Intrinsic & Extrinsic Age
Reported in almost every publication, this helps us most decide clinical recommendations

Telomere Length
More correlated and predictive of outcomes than regular telomere length but not clinically relevant yet.

Exclusively Licensed from Duke for All Verticals
Contains Immune Cell Subset Novel IP

The DunedinPACE is the Most Predictive Clock

Endpoint	DunedinPACE			PhenoAge			GrimAge			DunedinPACE		
	HR	95% CI	p-value	HR	95% CI	p-value	HR	95% CI	p-value	HR	95% CI	p-value
Overall PACE	1.24	[1.14, 1.34]	<0.0001	1.00	[-0.01, 0.01]	0.9999	1.00	[-0.01, 0.01]	0.9999	1.24	[1.14, 1.34]	<0.0001
PhenoAge	1.00	[-0.01, 0.01]	0.9999	1.00	[-0.01, 0.01]	0.9999	1.00	[-0.01, 0.01]	0.9999	1.00	[-0.01, 0.01]	0.9999
GrimAge	1.00	[-0.01, 0.01]	0.9999	1.00	[-0.01, 0.01]	0.9999	1.00	[-0.01, 0.01]	0.9999	1.00	[-0.01, 0.01]	0.9999

Key HRs: Overall PACE HR = 1.24 [1.14-1.34], CVD HR = 1.18 [1.05-1.34], Katz ADL IRR = 1.27 [1.02-1.58], Nagi ADL IRR = 1.26 [1.02-1.54], Rosow-Breslau ADL IRR = 1.27 [1.08-1.50]; associations with stroke were similar to unadjusted models (HR = 1.33 [1.05-1.69]). Results for all models are reported in [Supplementary File 16](#).

Thus, Dunedin PACE adds incremental prediction over and above all clocks studied here."

TruDiagnostic

The DunedinPACE is the Most Predictive Clock


Even just being slightly above an aging rate of 1 biological year/chronological year can increase your risk of death by 56% in the next 7 years and increase your risk of a chronic disease diagnosis by 54% over the next 7 years.



In fact, patients who were considered fast agers, were 16% more likely to die and 23% more likely to develop a chronic disease. That means they were 65% more likely to die in our cohorts than those at normal or slow aging.

TruDiagnostic


One-leg Balance Test





Pace of Aging 26 - 45 R = 0.36, p < 0.001


Each plotted point represents 20 study members


TruDiagnostic
















The DunedinPACE has the Highest Correlations to Quality of Life Metrics



TruDiagnostic

The DunedinPACE is the Most Precise Clock



TruDiagnostic


The DunedinPACE is modifiable by interventions which we already KNOW improve healthspan and lifespan

***This is EXTREMELY important

Biomarker Criteria	Horvath epigenetic age	Hannum epigenetic age	GrmAge	PhenoAge	DunedinPoAm
DNA Methylation Biomarker Calibrated to Detect:	Chronologic Age	Chronologic Age	Biomarkers, Smoking, Death	Phenotypic Age	Pace of Aging (change)
Feasible for use in a clinical trial in older adults?	✓	✓	✓	✓	✓
Robustly associated with chronological age across independent cohorts?	✓	✓	✓	✓	✓
Predicts age-related change in function, chronic disease, or death?	✓	✓	✓	✓	✓
Responsive to interventions that beneficially affect the biology of aging?	—	—	—	—	—

TruDiagnostic

The DunedinPACE is modifiable by interventions which we already KNOW improve healthspan and lifespan




Change from baseline to 12- and 24-month follow-up in DunedinPACE (Pace of Aging) measures of aging in ad libitum (AL) and caloric restriction (CR) groups in CALERIE Trial (Waziry, R., et al.)

TruDiagnostic

The DunedinPACE is modifiable by interventions which we already KNOW improve healthspan and lifespan


CALERIE RCT of caloric restriction (N=197)



TruDiagnostic

The DunedinPACE is modifiable by interventions which we already KNOW improve healthspan and lifespan

ADNI Study of ADRD (N=649)



TruDiagnostic

**Leukocyte Telomere Length:
An Epigenetic Predictor**
(accuracy/caveats)

Leukocyte DNAmTL (telomere length) is applicable across the entire age spectrum and is more strongly associated with age than measured leukocyte TL (LTL) ($r \sim -0.75$ for DNAmTL versus $r \sim -0.35$ for LTL).

DNA methylation-based estimator of telomere length

Alex T. Lu,¹ Anne Seeboth,² Pui-Chien Tsai,^{3,4,5} Dianliang Sun,^{6,7} Austin Quach,¹ Alex P. Reiner,⁸ Charles Kooperberg,⁹ Luigi Ferrucci,⁹ Lifano Hou,¹⁰ Andrea A. Baccarelli,¹¹ Yun Lu,¹² Sarah E. Harris,^{13,14} Anne Corley,^{15,16} Annie Santos,^{15,16} Ian J. Deary,^{15,16} James D. Stewart,¹⁷ Eric A. Whitson,^{15,16} Themistokles L. Assimes,^{17,18} Wei Chen,⁷ Shengou Li,¹⁹ Massimo Mangino,⁷ Jordana T. Bell,³ James G. Wilson,²⁰ Abraham Aviv,²¹ Riccardo E. Marioni,^{2,13} Kenneth Raj,^{22*} and Steve Horvath^{19,23*}

**Epigenetic Leukocyte Telomere Length:
A More Accurate Predictor of Health Outcomes**

Leukocyte DNAmTL outperforms LTL in predicting

- Time-to-death ($p=2.5E-20$)
- Time-to-coronary heart disease ($p=6.6E-5$)
- Time-to-congestive heart failure ($p=3.5E-6$)
- Association with smoking history ($p=1.21E-17$)


DNAmTL is not only an epigenetic biomarker of replicative history of cells, but a useful marker of age-related pathologies that are associated with it

Leukocyte Telomere Length:

Shorter LTL is associated with increased EEAA ($r = -0.16$, $p = 3.1 \times 10^{-6}$). LTL is inversely related to proportions of memory CD8+ T cells ($p = 4.04 \times 10^{-16}$) and positively related to proportions of naive CD8+ T cells.

Blood that contains more memory CD8+ T cells and less naive CD8+ T cells would display a relatively shorter LTL and older DNA methylation age.

EEAA is highly predictive of all-cause mortality. Epigenetic mortality risk is strongly associated with telomere length.



Better together!


From the previous slide we know that chronological clocks, phenotypic clocks, and telomere length don't really correlate well with each other meaning they represent different processes.

"Evidence that TL and epigenetic clock estimates are independent predictors of chronological age and mortality risk was obtained in the study by Marioni et al. (2018) performed in two Scottish cohorts aged from 70 to 90 years.

In both cohorts studied, combined whole-blood TL and DNAm age explained more variance in age than each of them individually. In a combined cohort analysis, TL and DNAm age explained 2.8 and 28.5% of the variance in age, respectively, and jointly they explained 29.5%. Also in a combined cohort, one standard deviation increase in a baseline DNAm age was associated with a 25% increased mortality risk ($p < 0.0001$) while in the same model, one standard deviation increase in a baseline TL was independently associated with an 11% reduced mortality risk only ($p = 0.05$."

TruDiagnostic

The Difficulty in Reversing Aging:



DOI: 10.1371/journal.pbio.1001047.g001
Figure 1. Physical and Actual Escape Velocities
 Remaining life expectancy follows a similar trajectory whether one walks off a cliff or reverses age; the slow walker differs, but most progressively worsens with time. Slight mitigation of this cohorter by jet propulsion or by rejuvenation therapies merely postpones the outcome, but sufficiently impedes deterioration outcomes. The degree of gene or protein and metabolite decreases the individual from a rocky end. Numbers denote plausible ages, at the time five-generation rejuvenation therapies arrive, of people following the respective trajectories.

At the moment, it is VERY difficult to make major differences in aging!

Aging is extremely complex and multi-factorial. While we would love to see these age markers decrease every time we test. This isn't likely. If it were, we would be immortal

We can SLOW this process for now.

TruDiagnostic

Case Studies



Epigenetic Age Reduction Case Study: Patient #1


Patient#1 first epigenetic age test was in September of 2020.

At a chronological age of 37.26, the patient had the following age metrics on the first test.

Patient implemented, Fisetin, Vitamin D, Caloric Restriction, and Stress maintenance exercises.



Patient Initials	Kil Number	Chronological Age	TL	PC Extrinsic	PC Intrinsic	DunedinPACE
ANON-1	#1	37.26	7.186274	36.65	51.746611	1.01




The DunedinPACE is modifiable by interventions which we already KNOW improve healthspan and lifespan

***This is EXTREMELY important

Biomarker Criteria	Horvath epigenetic age	Hannum epigenetic age	GrimAge	PhenoAge	DunedinPcAm
DNA Methylation Biomarker Calibrated to Detect:	Chronologic Age	Chronologic Age	Biomarkers, Smoking, Death	Phenotypic Age	Pace of Aging (change)
Feasible for use in a clinical trial in older adults?	✓	✓	✓	✓	✓
Robustly associated with chronological age across independent cohorts?	✓	✓	✓	✓	✓
Predict age-related change in function, chronic disease, or death?	✓	✓	✓	✓	✓
Responsive to interventions that beneficially affect the biology of aging?	---	---	---	---	---

TruDiagnostic

The DunedinPACE is modifiable by interventions which we already KNOW improve healthspan and lifespan



Change from baseline to 12- and 24-month follow-up in DunedinPACE (Pace of Aging) measures of aging in ad libitum (AL) and caloric restriction (CR) groups in CALERIE Trial (Waziry, R., et al.)

TruDiagnostic

Work-related stress and well-being in association with epigenetic age acceleration: A Northern Finland Birth Cohort 1966 Study

Major Lifestyle Recommendations: Reduce Stress

Our results indicate few significant associations between work stress indicators and epigenetic age acceleration, limited to a range of ± 2 years.

PhenoAgeAA was associated with job strain active work ($\beta=-1.301$, 95%CI -2.391, -0.212), slowing aging of less than 1.5 years, and working as white-collar slowed aging six months (CrimAgeAA $\beta=-0.683$, 95%CI -1.264, -0.102) when compared to blue collars.

Association was found for working for more than 40 hours per week that increased the aging over 1.5 years. (HorvathAA $\beta=2.058$, 95%CI 0.517, 3.599, HannumAA $\beta=1.567$, 95%CI 0.415, 2.719).

Rate of aging: DunedinPoAm
 For the unadjusted fit, we found a statistically significant association for being working as white-collar: $\beta=-0.015$, 95%CI -0.026, -0.003 compared with being blue-collar. Effort ($\beta=0.009$, 95%CI 0.001, 0.016) and evening shift increase the rate of aging of less than 2% ($\beta=0.015$, 95%CI 0.001, 0.03).

Exposure	Description
Job stress	Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job stress. High job stress was defined as a score of 20 or higher on the job stress scale of the Job Stress Scale (JSS) questionnaire (100 items, 0-100). Job stress was defined as a score of 20 or higher on the JSS questionnaire.
Job strain	Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job strain. High job strain was defined as a score of 20 or higher on the job strain scale of the JSS questionnaire (100 items, 0-100). Job strain was defined as a score of 20 or higher on the JSS questionnaire.
Job strain active work	Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job strain active work. High job strain active work was defined as a score of 20 or higher on the job strain active work scale of the JSS questionnaire (100 items, 0-100). Job strain active work was defined as a score of 20 or higher on the JSS questionnaire.
Job strain passive work	Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job strain passive work. High job strain passive work was defined as a score of 20 or higher on the job strain passive work scale of the JSS questionnaire (100 items, 0-100). Job strain passive work was defined as a score of 20 or higher on the JSS questionnaire.

Table 1. Job exposure definitions

Exposure

Description

Job stress
 Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job stress. High job stress was defined as a score of 20 or higher on the job stress scale of the Job Stress Scale (JSS) questionnaire (100 items, 0-100). Job stress was defined as a score of 20 or higher on the JSS questionnaire.

Job strain
 Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job strain. High job strain was defined as a score of 20 or higher on the job strain scale of the JSS questionnaire (100 items, 0-100). Job strain was defined as a score of 20 or higher on the JSS questionnaire.

Job strain active work
 Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job strain active work. High job strain active work was defined as a score of 20 or higher on the job strain active work scale of the JSS questionnaire (100 items, 0-100). Job strain active work was defined as a score of 20 or higher on the JSS questionnaire.

Job strain passive work
 Report an individual's perception of a psychosocial work environment that under a high level of job demands combined with a low level of control, resulting in job strain passive work. High job strain passive work was defined as a score of 20 or higher on the job strain passive work scale of the JSS questionnaire (100 items, 0-100). Job strain passive work was defined as a score of 20 or higher on the JSS questionnaire.

Phenotypic Age (PhenoAge)
 Derived from a linear regression model of the following variables: albuminuria, systolic blood pressure, waist circumference, waist-hip ratio, hemoglobin A1c, and total cholesterol.

Phenotypic Age (Hannum)
 An additive 50 CpG model, based on a regression model of the following variables: albuminuria, waist circumference, waist-hip ratio, hemoglobin A1c, and total cholesterol.

Phenotypic Age (Horvath)
 Phenotypic age is the sum of the following variables: hemoglobin A1c, waist circumference, waist-hip ratio, hemoglobin A1c, and total cholesterol.

Work-related variables
Work hours
 Working hours were categorized into three groups: less than 40 hours per week, 40-49 hours per week, and 50 or more hours per week.

Work intensity
 Work intensity was categorized into three groups: low, medium, and high.

Work engagement
 Work engagement was categorized into three groups: low, medium, and high.

Work effort
 Work effort was categorized into three groups: low, medium, and high.



Epigenetic Age Acceleration: Lifetime Stress

Cumulative stress was associated with accelerated GrimAge ($P=0.038$) and stress-related physiologic measures of adrenal sensitivity (Cortisol/ACTH ratio) and insulin resistance (HOMA). After controlling for demographic and behavioral factors, HOMA correlated with accelerated GrimAge ($P=0.016$).

Remarkably, psychological resilience factors of emotion regulation and self-control moderated these relationships. Emotion regulation moderated the association between stress and aging ($P=0.032$), such that with worse emotion regulation, there was greater stress-related age acceleration, while stronger emotion regulation prevented any significant effect of stress on GrimAge. Self-control moderated the relationship between stress and insulin resistance ($P=0.0073$), with high self-control blunting this relationship. In the final model, in those with poor emotion regulation, cumulative stress continued to predict additional GrimAge Acceleration even while accounting for demographic, physiologic, and behavioral covariates.

These results demonstrate that cumulative stress is associated with epigenetic aging in a healthy population, and these associations are modified by biobehavioral resilience factors.

Psychological and biological resilience modulates the effects of stress on epigenetic aging

Figure 3. Psychological resilience factors moderate the effects of cumulative stress on GrimAge Acceleration and physiologic stress pathways.

A GrimAge Acceleration (Years) vs. Cumulative Stress (Cohen's D). Interaction with Emotion Regulation: $P=0.032$. High emotion regulation blunts the relationship.

B HOMA vs. Cumulative Stress (Cohen's D). Interaction with Self-control: $P=0.0073$. High self-control blunts the relationship.

C Cortisol/ACTH ratio vs. Cumulative Stress (Cohen's D). Interaction with Emotion Regulation: $P=0.032$. High emotion regulation blunts the relationship.



Epigenetic Age Acceleration: Relaxing Treatment

Moreover, we compare DNAmAge with another mechanism of biological age, leucocyte telomere length (LTL) and telomerase. DNAmAge is reduced after training in healthy subjects ($p=0.053$), but not in patients. LTL is preserved after intervention in healthy subjects, while it continues to decrease in patients ($p=0.05$).

The conventional negative correlation between LTL and chronological age becomes positive after training in both patients ($p=0.01$) and healthy subjects ($p=0.05$), in our subjects, DNAmAge is not associated with LTL.

Our findings would suggest that intensive relaxing practices influence different aging molecular mechanisms, i.e., DNAmAge and LTL, with a rejuvenating effect.

Our study reveals that DNAmAge may represent an accurate tool to measure the effectiveness of lifestyle-based interventions in the prevention of age-related diseases.

Exploiting Epigenetic Age in Response to Intensive Relaxing Training: A Pilot Study to Slow Down Biological Age




Epigenetic Age Reduction Case Study: Patient #1

Patient#1 first epigenetic age test was in September of 2020.

At a chronological age of 37.26, the patient had the following age metrics on the first test.

Patient implemented, Fisetin, Vitamin D, Caloric Restriction, and Stress maintenance exercises.



TruDiagnostic

Epigenetic Age Reduction Case Study: Patient #1

Patient#1 first epigenetic age test was in September of 2020.

At a chronological age of 37.26, the patient had the following age metrics on the first test.

The patient tested again, 93 years later.

Patient implemented, Fisetin, Vitamin D, Caloric Restriction, and Stress maintenance exercises.

We saw Age reduction in each marker!

Patient Initials	Kit Number	Chronological Age	TL	PC Extrinsic	PC Intrinsic	DunedinPACE
ANON-1	#1	37.26	7.18	36.65	51.746611	1.01
ANON-2	#2	38.19	7.20	37.19	51.96	.95
Delta		.93	.02	0.54	0.213389	-0.06

TruDiagnostic

Diabetes has many markers in the Methylome

Loci Specific methylation can yield clinical insight.

> Epigenetics. 2019 Jul 2;1(7):482-8. doi: 10.1080/105802294.2019.1718418. Epub 2019 May 6.

DNA methylation of loci within ABCG1 and PHOSPHO1 in blood DNA is associated with future type 2 diabetes risk

Tarantini Dapini ¹, Trivandita Taneja ², X. H. ³, Dilar Altunbas ⁴, Alexander Perleberg ¹, Ben-Adar Shoshita ⁵, Vincenzo Di Stefano ⁶, Jussu Phipps ⁷, Alan Wang ⁸, Lutz Grunert ⁹, Emma Nilsson ¹⁰, Charlotte Ling ¹¹

> Epigenetics. 2019 Jun;1(10):885-897. doi: 10.2191/epi-2018-0168. Epub 2019 Jun 6.

Critical evaluation of the DNA-methylation markers ABCG1 and SREBF1 for Type 2 diabetes stratification

Christin Krause ¹, Helen Savari ², Cathleen Gellner ³, Martina Dreha ⁴, Alexander T. El-Gemmal ⁵, Stefan Walter ⁶, Olena Ocher ⁷, Falken Kilger ⁸, Ulrike M. Kilian ⁹, Meike Karsten ¹⁰, Christina Klein ¹¹, Georg E. Brunner ¹², Oliver Mann ¹³, Heidekl Lohmann ¹⁴, Heidekl Kretschmer ¹⁵

TruDiagnostic

Diabetes Risk Report

DNA methylation at the **ABCG1** locus **cg06500161** in blood DNA was associated with a 9% increased risk for future T2D (OR = 1.09, 95% CI = 1.02-1.16, P-value = 0.007, Q-value = 0.016)

DNA methylation at the **PHOSPHO1** locus **cg02650017** in blood DNA was associated with a decreased risk for future T2D (OR = 0.95, 95% CI = 0.75-0.95, P-value = 0.006, Q-value = 0.018) after adjustment for age, gender, fasting glucose, and family relation.

A DNA methylation score of 5.0% or greater at the **PHOSPHO1** locus **cg02650017** in blood DNA was associated with a 15% decreased risk for future type 2 diabetes occurrence.

Response to Caloric Restriction

CpG site	Gene	B-allele frequencies	Your score	Response Status
1. cg11500690	PKNO3	0.072	0.21	Hypomethylated
2. cg13131312	PKNO3	0.110	0.113	Hypomethylated
3. cg11433508	PKNO3	0.165	0.181	Hypomethylated
4. cg03013582	PKNO3	0.120	0.117	Hypomethylated
5. cg08081156	PKNO3	0.162	0.287	Hypomethylated
6. cg04082052	PKNO3	0.324	0.325	Hypomethylated
7. cg0447150	PKNO3	0.490	0.494	Hypomethylated
8. cg10229418	PKNO3	0.232	0.250	Hypomethylated
9. cg17159913	PKNO3	0.252	0.253	Hypomethylated
10. cg04790391	PKNO3	0.355	0.359	Hypomethylated
11. cg08481772	PKNO3	0.418	0.417	Hypomethylated

Possible Outcomes:

Non-Responder: CpG loci around your PKNO3 gene are generally under-methylated. Your response to caloric restriction alone is low, so it is not likely to be the most effective form of weight loss.

Intermediate Responder: CpG loci around your PKNO3 gene are in the normal range. Using a calorie deficit diet for weight loss may work, but it may not be as successful as other therapies.

Full Responder: CpG loci around your PKNO3 gene have been hypermethylated, so caloric restriction as a method of weight loss should be very effective.

YOUR EPITYPE:

Based on your methylation, you are a **non-responder** for weight-loss interventions using caloric restriction.

TruDiagnostic™

Smoking and Drinking Reports

This can be used as a good metric to measure patient adherence.

Smoking & Disease Risk

AHRH (cg00779021)
Average Methylation %

Your Epitope Value: 0.22

Score Range: 0.22 to 0.32

The report that follows outline measure risk for the epigenome is based on the level of methylation at the AHRH gene from epigenetics.

Your DNA methylation score was 0.22, which is in the low range. Having a low methylation score goes with the status of non-smokers. Getting your DNA methylation report to help you in developing smoking-related conditions.

Alcohol Consumption and DNA Methylation

Non-Smoker/No Alcohol | **Non-Smoker/Alcohol** | **Smoker/No Alcohol** | **Smoker/Alcohol**

Get your methylation score and compare your drinking status on 0-3. **Think For Yours!** With our custom methylation risk score, you are in the **low** percentile. This means your score is higher than 80% of the population and low risk.


There are several self-reported smoking or "No Smoker" score associated to have in drinking status and has been combined with this "No Smoker" status.

TruDiagnostic™

Inflammatory Marker Reports:


The utility of this proxy phenotype is particularly evident for studies in which DNAm is quantified but for which IL-6 is unavailable. Additionally, utilizing a composite methylation score that integrates information from multiple sites could potentially provide a more reliable estimate of chronic inflammation, allowing for clearer insight when assessing its relation to health outcomes.

DNAm is a proxy for the proinflammatory cytokine IL-6 that associates with pertinent health, lifestyle, and cognitive outcomes



DNA methylation signatures of C-reactive protein associations with structural neuroimaging measures and major depressive disorder

© Claire Greaves, Hany Sami, Amy J. Bastiaan, Shuang Li, C. Daniel Harlow, A. Harvey, Miriam C. Burke, Emma L. Heister, Mark J. Adams, Stephen M. Lawrie, Kathryn L. Best, © Karen M. Walker, Steven W. Morris, David J. Reardon, © Jason M. Wessely, Douglas Smeeth, Gordon D. Williams, Anna-Louise Smith, Andrew Campbell, Riccardo E. Polidori, Simon R. Cox, Jonathan Cavanagh, © Andrew P. McIntosh, Heather C. Whalley
doi: <https://doi.org/10.1101/2020.08.12.20175987>



DNA methylation signature of chronic low grade inflammation and its role in cardio-respiratory diseases

© Claire Greaves, Hany Sami, Amy J. Bastiaan, Shuang Li, C. Daniel Harlow, A. Harvey, Miriam C. Burke, Emma L. Heister, Mark J. Adams, Stephen M. Lawrie, Kathryn L. Best, © Karen M. Walker, Steven W. Morris, David J. Reardon, © Jason M. Wessely, Douglas Smeeth, Gordon D. Williams, Anna-Louise Smith, Andrew Campbell, Riccardo E. Polidori, Simon R. Cox, Jonathan Cavanagh, © Andrew P. McIntosh, Heather C. Whalley
doi: <https://doi.org/10.1101/2020.08.12.20175987>

Abstract
We performed a multi ethnic, Engageome Wide Association study on 22,776 individuals to assess the DNA methylation signature of chronic low grade inflammation as measured by C-reactive protein (CRP). We find a 2,320 independent differentially methylated loci associated with CRP. These CpG sites show consistent associations across ethnicities, and are primarily located in enhancers, depleted in CpG islands. These genomic loci are predominantly shared by transcription factor binding sites and genomic enhancer regions. Mendelian randomisation analysis suggests altered CpG methylation is a consequence of increased blood CRP levels. Mediation analysis reveals obesity and smoking as important underlying driving factors for changed CpG-methylation. Finally, we find that an activated CpG signature significantly increases the risk for cardiovascular diseases and COPD.
